

Safari Elephant B-S1 | 1-section table



 Box dimensions: **220 x 136 x 75 cm** Weight: **+90 kg**

 Length: **200 cm**

 Height: **50-99 cm**

 Width: **120 cm**

 Weight: **150 kg**

 Max load: **200 kg**

Features:

- 1-section Bobath and Voyta treatment table
- Extra wide surface
- Two-layer padding with non-flammable, bio-compatible and scratchproof upholstery in many colours
- Powder coated, extremely wide and stable steel frame
- Personal Authorization Safety System with 2 magnetic safety keys
- HalloTronic® System for electric high/low height adjustment from 50 to 99 cm with a capacity up to 200 kg
- Powder coated steel high/low height adjustment foot bar accessible from each side of the table
- Non-slip rubber feet with 1 cm height adjustment allowing for the levelling of the table
- F4 wheeling system - dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with non-slip rubber feet


”

Special mechanism of lifting the table on 4 castors is very fast and easy to operate and allows for full mobility. High functionality of tables allow us to perform various therapies. In my Rehabilitation Center we work with three tables manufactured by Meden-Inmed and they fulfill all our requirements.

Tomasz Wojdylak

Manual Therapist, Owner of Tomed Rehabilitation Center

Table articulations:



Ordering information:

B-S1.F0: 1-section massage and treatment table with basic wheeling system composed of 2 small wheels and 2 legs with non-slip rubber feet.

B-S1.F4: 1-section massage and treatment table with dual-lever retractable castor wheeling system composed of 4 directional castors with central brake and 4 legs with non-slip rubber feet.

